

Single-use Polyp Trap IFU

Not made with natural rubber latex

Intended use:

The single-use Polyp Trap Five chamber trap is used for suction retrieval of endoscopically removed polyps.



Warnings and Precautions:

- Consult the medical literature relative to techniques, complications and hazards prior to the performance of any endoscopic procedure.
- Endoscopic procedures should only be performed by persons having adequate training and familiarity with endoscopic techniques and polyp retrieval.
- Do not attempt to reuse, reprocess, refurbish, remanufacture or resterilize this device. Boke Endo Medical did not design this device nor is it intended to be reused, reprocessed, refurbished, remanufactured, or resterilized. Performing such activities on this disposable medical device presents a safety risk to patients (i.e. compromised device integrity, cross-contamination, infection).

Contraindications:

- Contraindications include, but are not limited to, those specific to any endoscopic procedure.

Prior to Use:



1. Prior to clinical use, inspect and familiarize yourself with the device. If there is evidence of damage (i.e. deformed, cracked vessel, missing component, or damaged packaging), **do not use this product and contact your local Product Specialist.**

Directions for Use (Specimen Retrieval):

1. Prior to clinical use, familiarize yourself with the device and read all the Instructions for Use.
2. Remove the Polyp Trap Five chamber trap from the package and confirm that the arrow on lid points to chamber 1 (See Figure 1).
3. Connect suction system tubing to the center port on the Polyp Trap five chamber trap lid (See Figure 1).
4. Connect Polyp Trap five chamber trap tubing to suction connector (barb) on endoscope (See Figure 2).
5. To retrieve resected polyp, apply suction via the endoscope suction valve until polyp is retrieved in the number 1 section of the polyp trap.
6. After the polyp is trapped, turn lid to the right (counterclockwise) so the arrow on the lid to points to the unnumbered open position on the jar (See Fig 2). This allows for continuous procedural suctioning.

7. To retrieve any additional polyps, the lid can be turned to chambers 2, 3 and 4, 5 respectively (****same patient only****).
8. To remove the Polyp Trap five chamber trap from the endoscope, disconnect the tubing from the endoscope suction connector (barb) and the suction system tubing.
9. Reattach the suction system tubing to the endoscope connector (barb).
10. Process specimens in accordance with institutional guidelines.

Note: This trap is not intended to be used as a transport container for specimens





Product Disposal:

After use, this product may be a potential biohazard which presents a risk of cross-contamination. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.